

# Cushion Covers

## External & Internal Incontinent Covers Standard with all Cushions



Rear View



Front View

### External Incontinent Cover

- Waterproof material surrounds cushion
- All-Way stretch waterproof material covers the cushion top
- Reverse lapped zipper ensures no fluid can pass through
- Labelled for easy identification
- Extra external covers are encouraged, and available through ADP

**NEW**

### Internal Incontinent Cover

- Waterproof material surrounds cushion
- All-Way stretch waterproof material covers the cushion top
- Designed to wick fluid away from the user
- Reverse zipper ensures fluid does not penetrate cushion
- Breathable material protects & maintains skin integrity
- Red for easy recognition
- Protects & maintains skin integrity
- Labelled for caregiver ease

#### Dry N'Cool Waterproof & Wicking Cover Option

- Designed to wick away up to 1 litre of fluid away from the user
- Breathable materials increases comfort
- Protects & maintains skin integrity

#### Wipe N'Dry Waterproof & Wipe-able Cover Option

- Designed to be wiped clean without removal of cover
- Breathable materials increases comfort
- Ideal for Long Term Care Facilities

**FRONT !**  
Please do not remove this cover  
**PUT ON THE EXTERNAL COVER !**