



# V.G. Arm Trough



## Composition

- Contoured trough padded in memory foam
- Includes foam hand pommel
- Slides securely onto armrest pad
- Adjusts to 31 fixed positions
- Single pivot point facilitates elevation
- Positioning strap secures arm comfortably
- Trough attached to polycarbonate tray

## Benefits & Applications

- Tray slides onto armrest for transfer ease
- Provides infinite arm-to-body angle adjustments
- Trough may be rotated toward the body to accommodate internal rotations
- Trough may be rotated away from the body to accommodate external rotations
- Trough may be mounted forward and/or back on tray to accommodate client's elbow angle
- Tray flips to side of wheelchair to allow for client transfer, without removal of trough/arm rest
- Hand pommel may be carved to meet user's individual tone
- Sturdy construction, will not warp
- No tools required for mounting and/or removal of tray
- Custom materials, sizes and shapes available