



## V.G. Half Lap Tray



### Composition

- Polycarbonate tray, 3/8" thick
- Removable top pad
- Washable cover
- Slides securely onto arm pad

### Benefits & Applications

- Allows for easy transfers
- Tray flips down flush with wheelchair wheels for transfer ease
- Tray may be used without pad, if required
- Sturdy construction, will not warp
- No tools required for mounting and/or removal
- Custom materials, sizes and shapes available



## V.G. Full Lap Tray



### Composition

- Polycarbonate tray, 1/4" or 3/8" thick
- Padded elbow region
- Straps securely to armrests
- Wraparound nylon belt secures behind the user

### Benefits & Applications

- Provides comfort for user with padded elbow region
- Sturdy construction, will not warp
- No tools required for mounting and/or removal
- Custom materials, sizes and shapes available