

# Special Edition Cushions

## Maintain Level Seating



### Composition


- Incorporates any Vision cushion and a Solid Seat Insert
- Seat Insert is adhered to the bottom of the selected cushion
- Contained by incontinent DryN'Cool or WipeN'Dry covers
- Composed of natural material that is resistant to warping

### Benefits & Applications

- Eliminates hammock effect caused by sling seat upholstery
- Provides increased positioning support
- Seat Insert cannot be removed/misplaced, as it's adhered to the cushion
- Maximizes the surface area available to the user
- Maintains a levelled seating surface
- Ensures the longevity of the product



# Drop Seats



### One Piece System

- Incorporates:
- The selected cushion
  - Attached universal insert
  - 4 drop clips

### Universal Fit

- Accommodate wheelchair cross braces
- Drop clips adjust to all wheelchair rails
- Ensures ease of use for the user and caregiver

### Drop Clips

- Available as 0", 1" and 2" drop clips
- Assist in achieving seat to floor heights
- Create seat angles to aid in positioning