

The **V.I.P.** Cushion

Pressure Relief Cushion



Composition

- Molded cushion ensures structural integrity and postural support
- Topped with *GELfoam* (gel infused soft visco memory foam) and high resiliency foam layers
- Pronounced contours
- Rounded front profile increases user comfort
- Includes 2 incontinent covers, including *DryN'Cool* cover:
 - Designed to wick fluid from user
 - Breathable materials increase comfort
 - Protects & maintains skin integrity
- Optional *WipeN'Dry* external cover available, easy to disinfect

Benefits

- Provides increased comfort & support with top layer of *GELfoam*
- Facilitates immersion to reduce pressure
- Prevents heat buildup
- Protects and maintains user's skin integrity
- Contours promote postural stability and correct positioning
- Maintenance free and caregiver friendly



Application

- Ideal for users that require moderate to high pressure redistribution
- Provides comfort and support and relief
- Appropriate for
 - Manual and Power wheelchair users
 - Extended use, 6-8 hours
 - Non-weight bearing users
- Available in standard and custom sizes to accommodate varying needs