

VerySOFT Deep Leg Trough Cushions

**Superb Pressure Relief
with Pronounced Profile**



2 Memory Foam Layers



Contour Shape

Composition

- Superior 2 layers of VerySOFT Memory foam Reduce Pressure & Increase Comfort
- High resiliency foam base cushion prevents bottoming out
- Deep Leg Troughs promote postural stability
- Rounded front 2¼" profile accentuate pummel and abductors
- Includes 2 incontinent covers


Benefits

- Protects & maintains user's skin integrity
- Deep leg troughs promote foot propulsion
- Maintains a 2¼" profile with pronounced pummel & abductors
- Maintenance Free Cushion
- Lightweight, will not impede movement
- Guaranteed not to bottom out for 30 months
- 3 Standard Levels of Softness: VerySOFT, VerySOFT'R, Very-VerySOFT, allow for customization of sensation without increase cost

Application


- Intended for clients that need moderate to high pressure redistribution
- Pronounced pummel and abductors aid in correcting and maintaining postural stability
- Low profile, 2¼", assists in achieving seat to floor heights
- Front profile promotes foot propulsion
- Ideal for clients that need a soft and supportive cushion
- Customizable in shape and size to accommodate varying needs, including leg length discrepancies, pelvic obliquities, etc.

Pressure Mapping of the VerySOFT Cushion



User on wheelchair sling upholstery

The red segment indicates a dangerous pressure area of 400mmHg



User sitting on VerySOFT cushion

The blue segments indicate safe pressure levels of 30mmHg