

VerySOFT Cushions

Superb Pressure Relief



2 Memory Foam Layers



Contour Shape

Composition

- Superior 2 layers of VerySOFT Memory foam Reduce Pressure & Increase Comfort
- High resiliency foam base cushion prevents bottoming out
- Gentle contours promote postural stability
- Rounded front profile
- Includes 2 incontinent covers


Benefits

- Protects & maintains user's skin integrity
- Maintenance Free Cushion
- Lightweight, will not impede movement
- Low 2 $\frac{3}{4}$ " profile
- 3 Standard Levels of Softness: VerySOFT, VerySOFT^R, Very-VerySOFT, allow for customization of sensation without increase cost
- Guaranteed not to bottom out for 30 months


Application

- For clients that need moderate to high pressure redistribution
- For clients that need a soft and supportive cushion
- Customizable in shape and size to accommodate varying needs, including leg length discrepancies, pelvic obliquities, etc.

Pressure Mapping of the VerySOFT Cushion



User on wheelchair sling upholstery
The red segment indicates a dangerous pressure area of 400mmHg



User sitting on VerySOFT cushion
The blue segments indicate safe pressure levels of 30mmHg