

Cushion Summary Sheet

- **GELfoam™**
 - VerySOFT gel infused visco elastic memory foam
 - Facilitates immersion to protect and maintain user's skin integrity
 - Significantly reduces pressure
 - Regulates temperature and prevents heat buildup
- Lightweight contoured cushions support self-propulsion and positioning
- Maintenance free and caregiver friendly cushions



VerySOFT™

- : Contoured cushion with pommel and abductors
- : Ideal for maintaining positioning and protecting skin integrity
- : Includes 2 layers of *GELfoam™* (Gel infused visco elastic memory foam)
- : 2 ¾" prior to compression
- : Appropriate for day long sitting (8+ hours/day)



VIP™

- : Contoured cushion with pommel and abductors
- : Top layer of *GELfoam™* (Gel infused visco elastic memory foam)
- : Includes layer of High Resiliency foam to aid in pressure absorption
- : 2½" prior to compression
- : Appropriate for long sitting (6-8 hours/day)



EPIC™

- : Contoured cushion with pommel and abductors
- : Ideal for foot propellers that require lower front profile
- : Includes top layer of *GELfoam™* (Gel infused visco elastic memory foam)
- : 2¼" prior to compression
- : Appropriate for long sitting (3-6 hours/day)



VITAL™

- : Contoured cushion with pommel and abductors
- : Top layer of Convoluted Foam
- : Promotes air circulation for breezy feel
- : Provides positioning support
- : 2 ¾" prior to compression
- : Appropriate for casual sitting (up to 3 hours/day)

